### Coffee (f221) IgE

Test ID:	2915	CPT:	86003
Clinical Significance			
Coffee (f221) IgE			





Methodology: Immunoassay-Analyte Specific Reagents

# ECH STICS AEORATORY MEDICINE

### Container Serum Separator Tube (SST®)

## Transport Temperature Room temperature

Specimen(s)	
0.3 mL serum	
Minimum Volume: 0.15 mL/	
allergen	

Specimen Stability			
Room temperature: 14 days			
Refrigerated: 14 days			
Frozen: 30 days			

Reject Criteria	

Days Performed	

#### **Collection Instructions**

\*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.