

## ANA Screen, with Reflex to Panel I

Test ID: 249

CPT: 86038\*

### Clinical Significance

ANA screen through a combination of IA & IFA methods

The ANA Reflex provides physicians with a cost-effective and medically justified approach to evaluating a patient with suspected rheumatologic disease. Eleven antibodies associated with specific rheumatologic disease entities are analyzed and resulted in sequential tiers until positive findings are reported. The test is not intended for the work-up of autoimmune hepatitis or other non-rheumatologic diseases

### Test Details

**Components:** If the ANA IA or IFA Screen is positive, it will reflex to Pattern & Titer (6002) and the following antibodies : dsDNA Ab , Anti-SSA Ab, Anti-SSB Ab, Smith Ab, Sm RNP Ab, Ribosomal P Ab, Scl-70 Ab, Centromere B Ab, Jo-1 Ab, RNP Ab, Chromatin Ab  
If ANA Screen, IFA is positive, then ANA Titer and Pattern will be performed at an additional charge (CPT code(s): 86039).  
if ANA Screen, IA is positive then antibodies will be performed at an additional charge:  
dsDNA (CPT code(s): 86225), Sm/RNP (CPT code(s): 86235), RNP (CPT code(s): 86235), Sm (CPT code(s): 86235), Chromatin (CPT code(s): 86235), SSA (CPT code(s): 86235), SSB (CPT code(s): 86235), Scl-70 (CPT code(s): 86235), Jo-1 (CPT code(s): 86235), Ribosomal P (CPT code(s): 83516) and Centromere B (CPT code(s): 86235).

### Methodology:

### Reference Range

Negative

### Container

Serum Separator Tube (SST)

### Transport Temperature

Refrigerated

### Specimen(s)

Serum

### Specimen Stability

Refrigerated: 7 days  
Frozen: 30 days

### Reject Criteria

Gross hemolysis • Grossly lipemic  
• Grossly icteric

### Days Performed

Mon-Sat  
\*Turnaround time 48 hours  
from time received

### Collection Instructions

\*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.