## Allergen, Blueberry (f288) IgE

| Test ID:   | 2568                    |                       | СРТ:                       | 86003 |
|--|-------------------------|-----------------------|----------------------------|-------|
| Clinical Significance  |                         |                       |                            |       |
| Common names:<br>Blueberry, European blueberry, Lowb         | oush blueberry, Highbus | h blueberry, Whinberr | ry, Whortleberry, Bilberry |       |
|  |                         |                       |                            |       |
|  |                         | L.T.D                 |                            |       |
| Test Details   |                         | Reference Rai         | 900                        |       |
| Components: Allergen, Blueberry (f28                         | Ω\ lgF                  |                       | nge<br>TM                  |       |
| Joinponents. Aneigen, Dideserry (120                         | AB                      | < 0.35 KU/L           | CH                         |       |
|  | IAGN                    | OST                   | ICS                        |       |
| AT TI<br><b>Methodology:</b> Immunoassay (IA)                | HE FOREFRONT D          | F LABORATORY          | MEDICINE                   |       |
|  |                         |                       |                            |       |
| Container  | Transport Ten           | nperature             | Specimen(s)                |       |
| Serum Separator Tube (SST)                                   | Refrigerated            |                       | Serum                      |       |
|  |                         |                       |                            |       |
| Specimen Stability   | <b>Reject Criteria</b>  |                       | <b>Days Performe</b>       | d     |
| Refrigerated: 7 days<br>Frozen: 30 days                      |                         |                       | Mon-Sat                    |       |
|  |                         |                       |                            |       |
| Collection Instructions                                      |                         |                       |                            |       |
| Separate serum from cells within 45 minutes of venipuncture. |                         |                       |                            |       |
| separate serum nom cens within 43 m                          | midles of verificative. |                       |                            |       |

\*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.