Allergen, Horse Dander (e3) IgE

7.mergen, morse Danaer (es, 192				
Test ID:	2603	СРТ	:	86003
Clinical Significance				
Latin name: Equus caballus Source material: Dander Family: Equidae Common names: Horse, Domestic Direct or indirect contact with ani		es sensitization.		
est Details		Reference Range		
Components:	AB	< 0.35 KU/L	H	
	IAGN		CS	
Methodology: Immunoassay (IA)	THE FOREFRONT OF	LABORATORY M	EDICINE	
	Tuesday Tours	wat wa	Cunnium au (a)	
Container	Transport Tempe	erature	Specimen(s)	
erum Separator Tube (SST)	Refrigerated		Serum	
Specimen Stability	Reject Criteria		Days Performed	
Refrigerated: 7 days Frozen: 30 days	ineject entend		Mon-Sat	
Collection Instructions				
Separate serum from cells within 4	5 minutes of venipuncture.			

*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.