Allergen,	Guinea	Piq E	pithelia	(e6	) IqE

Гest ID:	2606	СРТ:	86003
Clinical Significance			
Latin name: Cavia porcellus Source material: Epithelium and of Family: Caviidae Common names: Guinea pig, Cav Direct or indirect contact with an Animal allergens are major comp	y imal allergens frequently caus		
Seek Deteile		Defense Dance	
Components:	-AB	Reference Range < 0.35 KU/L	
AT  Methodology: Immunoassay (IA)	THE FOREFRONT OF	OSTICS LABORATORY MEDICINE	
Container Serum Separator Tube (SST)	Transport Tempe Refrigerated	Specimen(s) Serum	
Specimen Stability	Reject Criteria	Days Perform	ned
Refrigerated: 7 days Frozen: 30 days		Mon-Sat	
Collection Instructions			
Separate serum from cells within 4	15 minutes of venipuncture.		

\*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.