


Allergen, Apple (f49) IgE

Test ID: 2849

CPT: 86003

Clinical Significance

Distinguish between pollen related apple allergy and fruit allergy due to LTP sensitization
Improve the risk assessment using allergen components
Improve management of apple allergic patients.



Test Details

Components:

Methodology: Immunoassay (IA)

Reference Range

< 0.35 KU/LTM

Container

Serum Separator Tube (SST)

Transport Temperature

Refrigerated

Specimen(s)

Serum

Specimen Stability

Refrigerated: 7 days
Frozen: 30 days

Reject Criteria

Days Performed

Mon-Sat

Collection Instructions

Separate serum from cells within 45 minutes of venipuncture.

*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.