

## Allergen, Hazelnut, f17

Test ID: 3080

CPT: 86003

### Clinical Significance

Latin name: *Corylus avellana*

Source material: Shelled nuts

Family: Betulaceae (Corylaceae)

Common names: Hazel nut, Hazelnut, Filbert, Cobnut, Cob

The terms "Filbert" and "Hazel nut" are often used interchangeably for nuts from all plants in the genus *Corylus*, such as *C. silvestris*, *C. maxima* and *C. columna*.

### Test Details

#### Components:

**Methodology:** Immunoassay (IA)

### Reference Range

< 0.35 KU/L<sup>TM</sup>

### Container

Serum Separator Tube (SST)

### Transport Temperature

Refrigerated

### Specimen(s)

Serum

### Specimen Stability

Refrigerated: 7 days

Frozen: 30 days

### Reject Criteria

### Days Performed

Mon-Sat

### Collection Instructions

Separate serum from cells within 45 minutes of venipuncture.

\*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.