


## Allergen, Oyster, (f290) IgE

Test ID: 3508

CPT: 86003

### Clinical Significance

Latin name: Ostrea edulis  
Source material: Fresh whole oyster  
Family: Ostreidae  
Oyster and crustacea extracts show common antigenic structures.  
Allergen Exposure  
Oysters are eaten raw, cooked or smoked.  
Oyster sauce is used for flavouring of various dishes.  
Unexpected exposure  
Undeclared oyster sauce in various dishes.



### Test Details

#### Components:

**Methodology:** Immunoassay (IA)

### Reference Range

< 0.35 KU/L<sup>TM</sup>

### Container

Serum Separator Tube (SST)

### Transport Temperature

Refrigerated

### Specimen(s)

Serum

### Specimen Stability

Refrigerated: 7 days  
Frozen: 30 days

### Reject Criteria

### Days Performed

Mon-Sat

### Collection Instructions

Separate serum from cells within 45 minutes of venipuncture.

\*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.