## Allergen, Pine nut, (f253) IgE

Allergen, Pine nat, (j255) igt			
Test ID:	3512	СРТ:	86003
Clinical Significance			
Latin name: Pinus edulis Family: Pinaceae Common names: Pine nut, Pignole	es, Pinon nut, Pignola, Pinyor	nut, Pine kernals	
Components:	AB IAGN	Reference Range < 0.35 KU/L	
Methodology: Immunoassay (IA)	THE FOREFRONT OF	LABORATORY MEDICINE	
Container	Transport Temp	perature Specimen(s)	
Serum Separator Tube (SST)	Refrigerated	Serum	
Specimen Stability	Reject Criteria	Days Performe	ed
Refrigerated: 7 days Frozen: 30 days		Mon-Sat	
Collection Instructions			
	AF with large for a discount of		
Separate serum from cells within 4	is minutes of venipuncture.		

<sup>\*</sup>The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.