Allergen, Maple, Box-Elder (t1) IgE

Test ID:	3534	СРТ:	86003
Clinical Significance			
		e, Manitoba maple, Box elder maple, V ndo maple, Red river maple, Stinking as	
Test Details		Reference Range	
		ECH DSTICS AEORATORY MEDICINE	
Container	Transport Temper		(s)
Serum Separator Tube (SST)	Refrigerated	Serum	
Specimen Stability	Reject Criteria	Days Perf	ormed
Refrigerated: 7 days Frozen: 30 days		Mon-Sat	

Collection Instructions

Separate serum from cells within 45 minutes of venipuncture.

*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.