Allergen, House dust, Greer (h1) IgE

	3538	СРТ:	86003
Clinical Significance			
Latin name: Greer Labs., Inc Mite, cockroach and animal all	lergens are major componer	nts of house dust. For more information, s	ee Insects and Mites.
Test Details		Poforonco Pango	
Components:	LAB	Reference Range < 0.35 KU/L	1
Methodology: Immunoassay (IA	AT THE FOREFRONT A)	OF LABORATORY MEDICINE	
Container	Transport To	emperature Specim	en(s)
Serum Separator Tube (SST)	Refrigerated	Serum	
Specimen Stability	Reject Crite	ria Davs P	erformed
Refrigerated: 7 days Frozen: 30 days		Mon-Sa	
Collection Instructions			
Separate serum from cells withi	n 45 minutes of venipunctur	e.	

*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.