Allergy, Hymenoptera Panel

Test ID:	7913	CPT:	86003 x5
Clinical Significance			
C L.T.D			
	R		
Test Details		Reference Range	
	ee (i1), Paper wasp (i4), White-	See individual analytes	
faced hornet (i2), Yello	ow hornet (i5), Yellow jacket (i3)		
	DIAGN	OSTICS	
		LABORATORY MEDICINE	
Methodology: Immunoassay (IA)			
Container	Transport Temp	erature Specimen(s)	
Serum Separator Tube		Serum	
Specimen Stability Refrigerated: 7 days	Reject Criteria	Days Perform Mon-Sat	ied
Frozen: 30 days			
Collection Instructions			
Separate serum from cells within 45 minutes of venipuncture.			
'	·		
*The CPT codes prov	vided are hased on AMA guidance and are	for informational purposes only. CPT coding	is the
		ons regarding coding to the paver being billed	