

Allergy, Trees Panel, IgE

Test ID: 7973

CPT: 86003 X5

Clinical Significance



Test Details

Components: Maple (Box Elder), t1, Mountain Cedar, t6, Elm, t8, White Pine, t16, Melaleuca, t2

Methodology: Immunoassay (IA)

Reference Range

See individual analytes TM

Container

Serum Separator Tube (SST)

Transport Temperature

Refrigerated

Specimen(s)

Serum

Specimen Stability

Refrigerated: 7 days
Frozen: 30 days

Reject Criteria

Days Performed

Mon-Sat

Collection Instructions

Separate serum from cells within 45 minutes of venipuncture.

*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.