

Allergy, Mold Panel

Test ID: 92170

CPT: 86003 x 6

Clinical Significance



Test Details

Components: Alternaria tenuis (m6), Aspergillus fumigatus (m3), Cladosporium herbarum (m2), Mucor racemosus (m4), Penicillium notatum (m1), Stemphylium botryosum (m10)

Methodology: Immunoassay (IA)

Reference Range

See individual analytesTM

Container

Serum Separator Tube (SST)

Transport Temperature

Refrigerated

Specimen(s)

Serum

Specimen Stability

Refrigerated: 7 days
Frozen: 30 days

Reject Criteria

Days Performed

Mon-Sat

Collection Instructions

Separate serum from cells within 45 minutes of venipuncture.

*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.