# **Anti Nuclear Profile**

Test ID: 6004 CPT: 86225, 83516, 86235 X9

## **Clinical Significance**

The antinuclear antibody (ANA) test is commonly used to look for autoantibodies that attack components of cells nucleus, triggering autoimmune disorders like lupus, infections, and other autoimmune diseases.



#### **Test Details**

**Components:** dsDNA Ab , Anti-SSA Ab, Anti-SSB Ab, Smith Ab, Sm RNP Ab, Ribosomal P Ab, Scl-70 Ab, Centromere B Ab, Jo-1 Ab, RNP Ab, Chromatin Ab



Methodology: Immunoassay (IA)

## Reference Range

See individual analytes



Treenouslogy. Immunicussay (ii

#### Container

Serum Separator Tube (SST)

## Transport Temperature

Refrigerated

#### Specimen(s)

Serum

#### **Specimen Stability**

Refrigerated: 7 days Frozen: 30 days

#### **Reject Criteria**

Gross hemolysis • Grossly lipemic • Grossly icteric

### **Days Performed**

Mon-Sat

#### **Collection Instructions**

\*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.