

## ANA Screen w/Reflex to Panel II

Test ID: 251

CPT: 86038

### Clinical Significance

ANA screen through a combination of IA & IFA methods

### Test Details

**Components:** If ANA Screen, IFA is positive, then ANA Titer and Pattern will be performed at an additional charge (CPT code(s): 86039). If ANA Screen, IA is positive then antibodies will be performed at an additional charge: dsDNA (CPT code(s): 86225), RNP (CPT code(s): 86235), Sm (CPT code(s): 86235), Chromatin (CPT code(s): 86235), Scl-70 (CPT code(s): 86235), Jo-1 (CPT code(s): 86235), Ribosomal P (CPT code(s): 83516) and Centromere B (CPT code(s): 86235).

**Methodology:** Immunofluorescence, Immunoassay

### Reference Range

See Laboratory Report

### Container

Serum Separator Tube (SST)

### Transport Temperature

Room Temperature

### Specimen(s)

Serum

### Specimen Stability

Room temperature: 4 days  
Refrigerated: 7 days  
Frozen: 30 days

### Reject Criteria

Gross hemolysis, Grossly lipemic,  
Grossly icteric

### Days Performed

### Collection Instructions

\*The CPT codes provided are based on AMA guidance and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.