



Chlamydia trachomatis RNA, TMA

**Test Code:
 11361**

Transport temperature	Refrigerated or room temperature
Stability	Urine - 72 hours refrigerated; swab - 7 days room temperature
Turnaround time	24 - 72 hours
Accepted specimen	Urine container, Cobas® PCR swab, Orange Aptima® swab, White Aptima® swab, or ThinPrep® gynecological sample



Neisseria gonorrhoeae RNA, TMA

**Test Code:
 11362**

Transport temperature	Refrigerated or room temperature
Stability	Urine - 72 hours refrigerated; swab - 7 days room temperature
Turnaround time	24 - 72 hours
Accepted specimen	Urine container, Cobas® PCR swab, Orange Aptima® swab, White Aptima® swab, or ThinPrep® gynecological sample



Trichomonas vaginalis, PCR

**Test Code:
 90801**

Transport temperature	Refrigerated or room temperature
Stability	Urine - 72 hours refrigerated; swab - 7 days room temperature
Turnaround time	24 - 72 hours
Accepted specimen	Urine container, Orange Aptima® swab, White Aptima® swab, or ThinPrep® gynecological sample





Chlamydia/Neisseria gonorrhoeae RNA and Trichomonas

**Test Code:
 11364**

Transport temperature	Refrigerated or room temperature
Stability	Urine - 72 hours refrigerated; swab - 7 days room temperature
Turnaround time	24 - 72 hours
Accepted specimen	Urine container, Cobas® PCR swab, Orange Aptima® swab, White Aptima® swab, or ThinPrep® gynecological sample



Chlamydia/Neisseria gonorrhoeae RNA, PCR

**Test Code:
 11363**

Transport temperature	Refrigerated or room temperature
Stability	Urine - 72 hours refrigerated; swab - 7 days room temperature
Turnaround time	24 - 72 Hours
Accepted specimen	Urine container, Cobas® PCR swab, Orange Aptima® swab, White Aptima® swab, or ThinPrep® gynecological sample





LABTECHTM
DIAGNOSTICS
AT THE FOREFRONT OF LABORATORY MEDICINE

STD GUIDE

**CT/NG/Trich
Tests**

Chlamydia/Neisseria gonorrhoeae RNA, TMA, Throat

**Test Code:
70051**

Transport temperature	Refrigerated
Stability	30 days refrigerated
Turnaround time	7 - 10 business days
Accepted specimen	White Aptima [®] swab



Chlamydia/Neisseria gonorrhoeae RNA, TMA, Rectal

**Test Code:
16506**

Transport temperature	Refrigerated
Stability	30 days refrigerated
Turnaround time	7 - 10 business days
Accepted specimen	White Aptima [®] swab



SERVICES FOR EVERY DISCIPLINE



Herpes Simplex Virus 1&2-Specific Ab, IgG

**Test Code:
 164905**

Transport temperature	Refrigerated
Stability	14 days refrigerated
Turnaround time	24 hours
Accepted specimen	Serum separator tube (SST TM)



HIV 1/2 Antigen/Antibodies, with Reflexes

**Test Code:
 91431**

Transport temperature	Refrigerated
Stability	5 days refrigerated
Turnaround time	24 hours
Accepted specimen	Serum separator tube (SST TM)



Hepatitis Panel, Acute

**Test Code:
 322744**

Transport temperature	Refrigerated
Stability	14 days refrigerated
Turnaround time	24 hours
Accepted specimen	Serum separator tube (SST TM)



Syphilis Antibody Reflex RPR

**Test Code:
 90349**

Transport temperature	Refrigerated
Stability	7 days refrigerated
Turnaround time	24 - 48 hours
Accepted specimen	Serum separator tube (SST TM)

